

[WIP] Multi-Electrode Array Standard Operating Procedure



Hertie-Institut
für klinische Hirnforschung

by
Fabian Klopfer

Tübingen, 2022

Abstract

This standard operating procedure is supposed to provide information on the basic requirements and steps to be undertaken in order to record data from acute slices using the multi electrode array setup.

Much of the below is referring the reader to or adapted from the MultiChanel Systems MEA Application Note on Acute Hippocampus Slices from Rat or Mouse [8], the Axion Biosystems Acute Brain Slice Protocol [1], peer-reviewed protocols [4],[2, chapter 6] and the LAS interactive course on animal handling[11]. All manuals are linked in the citation, i.e. if you click the link, you'll download the document. The peer-reviewed protocols are very detailed, the notes by MCS however are also very detailed. Another source for this is the instructions that I received by Nikolas Layer and Daniela Miely.

Contents

1	Setup	3
2	Preparation	3
2.1	Artificial Cerebrospinal Fluid	3
2.2	Acute Slices	5
3	Recording Steps	5
4	Stimulation	5
5	Analysis	5

1 Setup

General documentation on this can be found in the MEA2100-System: Installation Guide brochure by MCS [10], CMOS-MEA5000-System: Quickstart [6], or in the respective user manuals [7, 5].

The Setup for the CMOS and the 256 MEA are very similar. The only difference is that the chip and the head stage have to be swapped (2100-256 and 5000-CMOS), as well as the recording software used.

Detailed handling instructions for the chips, like coating, cleaning and storage can be found in the MCS MEA manual [9].

Components:

- MEA-Chip (256MEA or CMOS-MEA; SCB or SCA for slices, CC for cell culture. For more details see the respective chip datasheets)
- MEA-Headstage (2100 for 256 or 5000 for CMOS)
- Interface board MCS-IFB-C 3.0 Multiboot
- MEA-to-interface board Connector cable (iX-industrial cable, type B)
- Temperature controller (TC02)
- Perfusion Cannula (PH01)
- Peristaltic perfusion system (PPS2)
- Video-Microscope
- Light source for the microscope
- Controls for the microscope

For chemical stimulation, an additional stand and a puffing electrode is necessary [3]. For measuring concentrations and DC two micromanipulators and electrodes are necessary, one of which must be prepared specifically.

TODO add pictures

2 Preparation

2.1 Artificial Cerebrospinal Fluid

Recipes for aCSF can also be found in [4, p. 10], [2, p. 122] and [8]. Example for recording FHM3 CSDs:

- 3,5 mM [K⁺] aCSF (recording aCSF): STOCK solution (10x)

Chemical	Quantity [g]
NaCl	72,466
KCl (for 35mM)	2,6093
MgCl ₂ * 6H ₂ O	2,03
NaH ₂ PO ₄	1,45
NaHCO ₃	21,85

Preparation of aCSF:

- 100ml of stock solution
- add 3,6g of Glucose (for 20mM)
- make up to 800ml volume with Millipore water and bubble with carbogen for at least 10min
- add CaCl₂ stock solution (2ml for recording)
- make up to a final volume of 1l
- final pH os 7.4 and osmolarity is 305mOSM

- 3,5mM [K⁺] aCSF (recording aCSF)

Chemical	Quantity
NaCl	7,2466g
KCl (for 3,5mM)	3,5ml (1M stock solution)
MgCl ₂ * 6H ₂ O	1ml (1M stock solution)
NaH ₂ PO ₄	0,145g
NaHCO ₃	2,185g
D-Glucose	3,6g
CaCl ₂	2ml (1M stock solution)

2.2 Acute Slices

For getting tissue, you should talk to your PI, if you don't have that already.

If you have an animal, the steps to produce slices are described in detail in [4, p. 10], [2, p. 122] and [8].

3 Recording Steps

More detailed instructions and protocols can be found in [4, p. 10], [2, p. 122], [1] and [8].

1. Turn on Devices: Pump, interface board, light source, video microscope, computer, controller for microscope/micromanipulators
2. Open Software: Pump, Temperature, Video Microscope, Multi Channel Experimenter
3. Set-up pump: connect the tubes as outlined in the Installation guide (printer version in one of the two metal cases below the table to the right of the setup). Test the perfusion system and make sure that the waste's input is lower than the perfusions output in the chamber!
4. Transfer and place slice: [8] mentions, that the slice should be placed on the electrodes without liquid for a better establishment of contact between the electrodes and the tissue. I.e. place it with the transfer pipette, soak up the extra liquid with paper tissue and quickly start the perfusion to avoid damage to the tissue.
5. Measure: In the MC Experimenter, draw the mea headstage into the main area. Drag the recording tool to the main area. Connect the digital data from the headstage to the recorder. Double click the recorder. In the bottom of the window there is a tab "data acquisition" click it to see the windows for the graphs to be recorded. Press start DAQ to get graphs. Press record to actually record data.
6. Clean & Power off: Clean the tubes with distilled water and hang them up such that they can dry out and don't mold. Transfer your data to an external hard drive. Turn off all devices that you turned on. Rinse the MEA with VE water and dry it, in case you haven't coated anything. Please look up the cleaning procedure necessary in the MEA manual [9].

4 Stimulation

TODO

1. Chemical
2. Electrical

5 Analysis

TODO

References

- [1] *Axion Bio - Culture Protocol - Acute Brain Slice*. Feb. 10, 2022. URL: <https://www.axionbiosystems.com/sites/default/files/resources/Axion%20Bio%20-%20Culture%20Protocol%20-%20Acute%20Brain%20Slice.pdf> (visited on 07/11/2023).
- [2] Mark Dallas and Damian Bell. *Patch clamp electrophysiology*. Springer, 2021.
- [3] Shahab Dehbandi et al. “Cortical spreading depression modulates synaptic transmission of the rat lateral amygdala”. In: *European Journal of Neuroscience* 27.8 (2008), pp. 2057–2065.
- [4] Kevin M Manz et al. “Patch-clamp and multi-electrode array electrophysiological analysis in acute mouse brain slices”. In: *STAR protocols* 2.2 (2021), p. 100442.
- [5] *MCS CMOS-MEA5000-System Manual.pdf*. Feb. 17, 2023. URL: https://www.multichannelsystems.com/sites/multichannelsystems.com/files/documents/manuals/MCS_CMOS-MEA5000-System_Manual.pdf (visited on 07/14/2023).
- [6] *MCS CMOS-MEA5000-System Quickstart Manual.pdf*. Feb. 10, 2022. URL: https://www.multichannelsystems.com/sites/multichannelsystems.com/files/documents/manuals/MCS_CMOS-MEA5000-System_Quickstart_Manual.pdf (visited on 07/14/2023).
- [7] *MCS MEA2100-System Manual.pdf*. Nov. 30, 2022. URL: https://www.multichannelsystems.com/sites/multichannelsystems.com/files/documents/manuals/MCS_MEA2100-System_Manual.pdf (visited on 07/14/2023).
- [8] *MEA Application Note on Acute Hippocampus Slice from Rat or Mouse*. Feb. 10, 2022. URL: https://www.multichannelsystems.com/sites/multichannelsystems.com/files/documents/applications/MEA%20Application%20Note_Acute%20Hippocampus%20Slice%20from%20Rat%20or%20Mouse.pdf (visited on 07/11/2023).
- [9] *MEA-System - MCS MEA Manual.pdf*. Jan. 23, 2023. URL: https://www.multichannelsystems.com/sites/multichannelsystems.com/files/documents/manuals/MCS_MEA_Manual.pdf (visited on 07/14/2023).
- [10] *MEA2100-System Installation Guide.pdf*. Feb. 10, 2022. URL: https://www.multichannelsystems.com/sites/multichannelsystems.com/files/documents/manuals/MEA2100-System_InstallationGuide.pdf (visited on 07/14/2023).
- [11] *Über LAS interactive*. July 11, 2023. URL: <https://www.las-interactive.de/index.php?id=6005#section0> (visited on 07/11/2023).